
## INSTRUCTIONS FOR USE

ETHICON

# SECURESTRAP™

Absorbable Strap Fixation Device

ETHICON SECURESTRAP™ 可吸收钉紧固装置

ETHICON SECURESTRAP™ 可吸收带固定装置

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### **ETHICON SECURESTRAP™ Absorbable Strap Fixation Device**

#### **Single Use Instrument**

Read all Directions, Contraindications, Warnings and Precautions prior to use.

This instructions for use provides directions for the proper use of the ETHICON SECURESTRAP™ Absorbable Strap Fixation Device.

This is not intended to be a comprehensive instruction manual for the performance of laparoscopic surgical procedures. Prior to use of this device, surgeons should be familiar with laparoscopic surgery and open hernia repair techniques.

#### **CAUTION**

Federal law (USA) restricts this device to sale by or on the order of a physician.

#### **DEVICE DESCRIPTION**

The ETHICON SECURESTRAP™ Absorbable Fixation Device contains either 12 or 25 synthetic absorbable straps, preloaded into the device's 36 cm long shaft. The instrument is designed for introduction and use through a 5 mm or larger laparoscopic port sleeve. Larger diameter sleeves will require use of a converter.

The ETHICON SECURESTRAP™ Absorbable Fixation Device straps are made of a blend of polydioxanone dyed with D&C Violet No. 2 and an L(-)-lactide/glycolide copolymer. The inserted length of the strap is 6.7 mm.

#### **ACTIONS**

In vitro studies indicate that the absorption profile of the polydioxanone and L(-)-lactide/glycolide copolymer is minimal in the first two weeks post implantation. The mechanism by which the polymer degrades is hydrolysis. After an initial reduction in the molecular weight of the strap, hydrolysis continues and a reduction in physical properties of the strap occurs. Significant absorption is seen after the critical tissue integration period has occurred, and the polymer is considered essentially absorbed at approximately 12-18 months.

#### **INDICATIONS FOR USE**

The ETHICON SECURESTRAP™ Absorbable Fixation Device is intended for fixation of prosthetic material to soft tissues in various minimally invasive and open surgical procedures, such as hernia repairs.

#### **CONTRAINDICATIONS**

- The device is not intended for use when prosthetic material fixation is contraindicated.
- Do not use the system on tissue that cannot be inspected visually for hemostasis.
- A minimum tissue thickness of 6.7 mm is required when applying the fastener over underlying bone, vessels, or viscera. If the total distance from the surface of the tissue to the underlying structure is less than the minimum tissue thickness, or may be compressed to a total distance less than the minimum tissue thickness, use of the device is contraindicated.
- This device should not be used in tissues that have a direct anatomic relationship to major vascular structures. This would include the deployment of fasteners in the diaphragm in the vicinity of the pericardium, aorta, or inferior vena cava during diaphragmatic hernia repair.

## WARNINGS


- This device is designed for single patient use only. Reuse or reprocessing of this device will lead to unpredictable strap performance and may lead to its failure and subsequent patient injury.
- This device is provided STERILE and is intended for use in a single patient. DISCARD AFTER USE. Do not reuse, reprocess or resterilize this device.
- The total distance from the surface of the tissue to the underlying bone, vessels, or viscera should be evaluated prior to application and should be a minimum of 6.7 mm.

## PRECAUTIONS

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

- Fixation method for any prosthetic device should be determined on the basis of accepted surgical technique, procedural requirements, and the instructions for use of the prosthetic material.
- Only physicians having adequate training and familiarity with endoscopic techniques should perform endoscopic procedures. A thorough understanding of the operating principles, risks versus benefits, and the hazards involved in utilizing an endoscopic approach is necessary to avoid possible injury to the user and/or patient.
- Verify mechanical and electrical compatibility of devices from different manufacturers prior to using them together in a procedure.
- Inspect fixation site after application to ensure hemostasis.
- Appropriate force should be applied to the device for application of straps. Excessive force may result in damage to the tissue, the material being fixated, or the device.


Note: Care should be taken in the application of external counter-pressure to avoid tissue compression between the end of the device and the location of counter-pressure of less than 6.7 mm.

## ADVERSE REACTIONS

Adverse reactions associated with the use of this device include transitory local irritation at the implant site and a transitory inflammatory foreign body response. As with all foreign devices, potential adverse events may also include inflammation and potentiation of procedure-related infection.

## SCHEMATIC VIEW

## INSTRUCTIONS FOR USE

1. Insert the ETHICON SECURESTRAP™ Absorbable Fixation Device through the appropriately sized laparoscopic port sleeve or larger port sleeve with the use of a converter (if necessary).
2. Place the tip of the device on the tissue or mesh, directly over the target site. The tip of the device may be used to facilitate positioning of the mesh.
3. While applying appropriate external counter-pressure, squeeze the handle completely, and then release. The application of counter-pressure helps ensure secure strap placement.
4. Repeat the above steps, 2-3, for placement of additional straps all along the length of the area to be fastened.
5. NOTE: If necessary, the straps may be dislodged from the mesh by cutting the strap with a laparoscopic scissor. Pull the mesh away from the strap and leave the strap in the tissue to resorb.
6. Following the deployment of the 12th or 25th strap, the handle will lock in the fully depressed position and the window on the device handle will be fully colored. The device should be removed and discarded at this time. If additional fixation points are required, obtain an additional device.
7. Upon completion of the application, remove the ETHICON SECURESTRAP™ Absorbable Fixation Device from the laparoscopic port sleeve.
8. The device may also be used in open surgical procedures.

## STORAGE


Store at room temperature. Avoid prolonged exposure to elevated temperatures. Do not expose to temperatures above 122°F (50°C).

## HOW SUPPLIED

The ETHICON SECURESTRAP™ Absorbable Fixation Device is available in single packages as sterile devices. The product is available in 6-up configurations.

Dispose of the device and packaging according to your facility's policies and procedures concerning biohazardous materials and waste.

## MEANING OF THE SYMBOLS

|   |   |   |  |
|---|---|---|--|
|   | Do Not Reuse  |   | Batch number   |
|  | Do Not Resterilize  |  | Upper Limit Temperature  |
|  | Do Not Use if Package is Damaged  |  | Catalogue number   |
|  | Use by – year & month   |  | Manufacturer   |
|  | Method of Sterilization – Ethylene Oxide  |  | CE-mark and Identification number of Notified Body. The product meets the essential requirements of the Medical Device Directive 93/42/EEC |
|  | See Instructions for Use  |  |  |
|  | CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician |   |  |
|  | Authorized Representative in the European Community   |   |  |

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© Ethicon, Inc. 2010
Johnson & Johnson Medical GmbH
ETHICON, LLC