ETHICON Harmonic®

Harmonic® HD 1000i Shears
Ciseaux Harmonic® HD 1000i
Harmonic® HD 1000i Schere
Forbici Harmonic® HD 1000i
Tesoura Harmonic® HD 1000i
Tijeras HD 1000i Harmonic®
Harmonic® HD 1000i-schaar
Harmonic® HD 1000i saks
Harmonic® HD 1000i -sakset
Ψαλίδι Harmonic® HD 1000i
Harmonic® HD 1000i-sax
Nožyce Harmonic® HD 1000i
Harmonic® HD 1000i olló
Nůžky Harmonic® HD 1000i
Nožnice Harmonic® HD 1000i
Harmonic® HD 1000i saks
Harmonic HD 1000i Makas
Ножницы Harmonic® HD 1000i
Foarfeca Harmonic® HD 1000i
Gunting HD 1000i Harmonic®
Dụng cụ cắt Harmonic® HD 1000i
Käärid Harmonic® HD 1000i
Harmonic® HD 1000i škéres
Harmonic® HD 1000i žirklés
Harmonic® HD 1000i ножици
Škarje Harmonic® HD 1000i
Škarje Harmonic® HD 1000i
Ножици Harmonic® HD 1000i
Маказе Harmonic® HD 1000i
Маказе Harmonic® HD 1000i
Harmonic® HD 1000i 手术剪

Please read all information carefully.








Failure to properly follow the instructions may lead to serious surgical consequences.

Important: This package insert is designed to provide instructions for use of the HARMONIC® HD 1000i Shears. It is not a reference to surgical techniques.

ETHICON

PART OF THE **ETHICON-ETHICON** FAMILY OF COMPANIES

Instructions, Instruccions, Gebrauchsanweisung, Istruzioni, Instruções, Instrucciones, Instructies, Brugervejledning, Ohjeet, Oðilycis, Anvisningar, Instrukcje, Utasítások, Návod k použití, Návod na použitie, Instruksjoner, Talimatlar, Инструкции, Instrucțiuni, Instruksi, Hướng dẫn, Kasutusjuhend, Instrukcijas, Instrukcija, Инструкции, Упуте, Navodila, Упутства, Упутства, Упутства, 使用说明

Indications

The HARMONIC HD 1000i Shears instrument is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The instruments can be used as an adjunct to or substitute for electrosurgery, lasers and steel scalpels in general, plastic, pediatric, gynecologic, urologic, thoracic, exposure to orthopedic structures (such as spine and joint space), sealing and transection of lymphatic vessels, and other open and endoscopic procedures. The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Energy button with Advanced Hemostasis.

Contraindications

- The instruments are not indicated for incising bone.
- The instruments are not intended for contraceptive tubal occlusion.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with ultrasonic devices include the potential for bleeding, tissue injury via mechanical or thermal damage, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, electrical shock, foreign body or magnetic resonance incompatibility, and property or environmental damage. Also, unintended harm, extended surgery, or altered surgical approach may result from issues related to device activation, damaged devices, electromagnetic interference, audible noise due to misassembly, misuse of the torque wrench, or an attempt to alter the device.

Device Description

The HARMONIC HD 1000i Shears instrument is a sterile, single patient use instrument used for dissection, grasping, coagulation, and cutting between the blade and clamp arm. It consists of an ergonomic handle with an integrated hand piece and two energy delivery buttons:

- 1) Energy button  - user can adjust power levels from 1-5.
- 2) Energy button with Advanced Hemostasis  - for large vessel sealing; user cannot adjust.

The instrument is available in two shaft lengths -- 20 cm and 36 cm.

An integrated audible and tactile mechanism in the handle indicates full trigger closure. The instrument has a clamp arm and coated curved blade that are designed to work through a 5 mm trocar, through a 5 mm reducer cap in a larger diameter trocar, or through an incision without the use of a trocar. The instrument shaft can be rotated continuously to facilitate visualization and access to targeted tissue. The two dashes on the instrument are intended to represent relative vessel size. The Energy button is indicated for vessels up to 5 mm in diameter. When the Energy button is used, cutting speed is the fastest. The Energy button with Advanced Hemostasis is designed for larger vessels and is indicated for vessels up to 7 mm in diameter. When the Energy button with Advanced Hemostasis is used, cutting speed is reduced and hemostasis is maximized. The instrument utilizes Adaptive Tissue Technology. This provides the generator with the ability to identify and monitor the instrument during use, which enables the generator to modulate and adjust its power output as well as provide audible feedback to the user as appropriate.

The HARMONIC HD 1000i Shears instrument is designed for use exclusively with the Generator G11 (GEN11) software version 2016-1 or later. Software revision can be found under "System Information" in the Generator G11 (GEN11) "Settings" menu. Refer to the Generator G11 (GEN11) Operator's Manual for more information.

Illustration and Nomenclature (Illustration 1)

- | | |
|-----------------------------|---|
| 1. Coated Blade | 5. Energy Button |
| 2. Clamp Arm and Tissue Pad | 6. Energy Button with Advanced Hemostasis |
| 3. Shaft | 7. Handle |
| 4. Rotation Knob | 8. Trigger |

Transport and Storage Conditions

Temperature: -22°C to +60°C
Relative Humidity: 10-80%

Warning and Caution Statements


WARNING: A Warning statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

Caution: A Caution statement alerts the user of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property. It may also be used to alert against unsafe practices. This includes the special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Instructions for Use

Verify compatibility of all instruments and accessories prior to using this instrument (refer to **Warnings and Precautions**).

- 1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
- 2 Connect the instrument to the generator and turn the generator power on.
- 3 Select the desired power level using the INCREASE and DECREASE buttons on the generator touchscreen. Only the power level for the Energy button can be adjusted (1-5). The system defaults to power level 5. (Illustration 2) For greater tissue cutting speed use a higher generator power level, and for greater coagulation use a lower generator power level. The amount of energy delivered to the tissue and resultant tissue effects are a function of many factors, including the power level selected, blade characteristics, grip force, tissue tension, tissue type, pathology, and surgical technique.
- 4 Close the clamp arm by closing the trigger, and insert the shaft through a trocar or incision (Illustration 3).
- 5 Position the tissue within the jaws at the desired location. The instrument's shaft can be rotated continuously using the rotation knob to facilitate visualization and access to target tissue.
- 6 Squeeze the trigger until it stops against the plastic handle (and a click is heard) to clamp targeted tissue between the jaws.
 - To achieve complete sealing, the trigger should be fully closed and the vessel fully contained between clamp arm and blade of device. An audible and tactile "click" indicates full trigger closure. To achieve full closure of the jaws of the device, squeeze the plastic trigger until you feel it stop against the plastic handle (plastic to plastic). Grip force needs to be maintained throughout the transection to keep trigger closed.
- 7 To activate the instrument blade, press one of the foot pedal switches or one of the energy buttons on the instrument.
 - Pressing the left foot pedal of the footswitch or the Energy button with Advanced Hemostasis on the instrument activates Advanced Hemostasis. When using the Energy button with Advanced Hemostasis, energy is not delivered unless the jaws are completely closed. This button activates an algorithm in the generator that, in conjunction with full trigger closure, allows for sealing larger vessels (up to 7 mm in diameter).
 - Pressing the right foot pedal of the footswitch or Energy button on the instrument activates the selected power level (1-5). The Energy button allows for sealing vessels up to 5 mm in diameter with full trigger closure and can enable other soft tissue applications (backcutting, scoring, drilling/otomy creation, etc.), where full trigger closure is not needed.

Foot Pedal	Button
Right	Energy Button 
Left	Energy Button with Advanced Hemostasis 

• The generator emits one of the audible tones listed in the table below to indicate when the instrument blade is first activated.

Tone	Button	Action
Repeating single tone	Energy Button	Generator On: Device is active
3 repeating, ascending tones	Energy Button with Advanced Hemostasis	Generator On: Device is active and in Advanced Hemostasis mode

- The generator changes to a second audible tone as Adaptive Tissue Technology regulates the delivery of energy. Thermal influences such as fluids or minimal to no tissue in the jaws may affect the presence or timing of the tone change. The tone change does not provide confirmation of tissue effect. When the second tone is heard, the situation should be assessed and the intended surgical action completed, such as gradual application of tension to facilitate transection. The secondary audible tone change is not a substitute for surgical experience.

Tone	Button	Action
High pitched, repeating single tone	Applies to both buttons	Adaptive Tissue Technology is active

WARNING: Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated (Illustration 4). Contact with staples, clips or other instruments while the instrument is activated may result in scratches on the blade, cracked or broken blades and premature blade failure.

WARNING: Do not use the Energy button with Advanced Hemostasis when energy application is desired prior to full closure of the jaws. Energy is not delivered using the Energy button with Advanced Hemostasis until the jaws are completely closed. Using the Energy button with Advanced Hemostasis without full trigger closure may result in lack of hemostasis.

WARNING: If activation is unintentionally stopped while sealing, maintain jaw closure and reactivate. Releasing the trigger while sealing may result in lack of hemostasis.

Caution: Keep the clamp arm open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm and distal shaft temperatures (Illustration 5).

WARNING: During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.

For optimal performance and to avoid tissue sticking, clean the instrument blade, clamp arm, and distal end of the shaft throughout the procedure by activating the instrument tip in saline. (Illustration 6)

WARNING: Do not clean the blade tip with abrasives. It can be wiped with a moist gauze sponge to remove tissue, if necessary. If tissue is still visible in the clamp arm, use hemostats to remove residue, taking care not to actuate the instrument (Illustration 7). Do not touch the instrument to the hemostats while activated. Scratches on the blade may lead to cracked or broken blades and premature blade failure.

- 8 Close the clamp arm by closing the trigger and remove the shaft from the trocar or incision.
- 9 Unplug the instrument from the generator.
- 10 Turn the generator **OFF** at the power switch.
- 11 Dispose of the instrument and cable in an appropriate container. No disassembly is required.

Warnings and Precautions

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.

- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse instruments in liquid unless the instruments are designed and labeled to be immersed.
- Verify compatibility with generators. Use device only with Ethicon Endo-Surgery Generator G11 (GEN11) software version 2016-1 or later. Software revision can be found under “System Information” in the Generator G11 (GEN11) “Settings” menu. Refer to the Generator G11 (GEN11) Operator’s Manual for more information.
- Avoid contact with any and all metal or plastic instruments or objects when the instrument is activated. Contact with staples, clips or other instruments while the instrument is activated may result in scratches on the blade, cracked or broken blades and premature blade failure.
- Do not use the Energy button with Advanced Hemostasis when energy application is desired prior to full closure of the jaws. Energy is not delivered using the Energy button with Advanced Hemostasis until the jaws are completely closed. Using the Energy button with Advanced Hemostasis without full trigger closure may result in lack of hemostasis.
- If activation is unintentionally stopped while sealing, maintain jaw closure and reactivate. Releasing the trigger while sealing may result in lack of hemostasis.
- Keep the clamp arm open when backcutting or while the blade is active without tissue between the blade and tissue pad to avoid damage to the tissue pad and increased blade, clamp arm and distal shaft temperatures.
- During benchtop testing of vessels >5 mm, the strongest vessel seals were achieved by allowing the Advanced Hemostasis mode to completely transect the targeted vessel.
- Do not clean the blade tip with abrasives. It can be wiped with a moist gauze sponge to remove tissue, if necessary. If tissue is still visible in the clamp arm, use hemostats to remove residue, taking care not to actuate the instrument. Do not touch the instrument to the hemostats while activated. Scratches on the blade may lead to cracked or broken blades and premature blade failure.
- In case of system failure, ensure the availability of the appropriate back up equipment relevant to the specific procedure.
- Audible high-pitched ringing, resonating from the blade, is an abnormal condition and an indicator that the blade is not operating properly. This may result in abnormally high shaft temperatures and user or patient injury.
- The instruments allow for the coagulation of vessels up to and including 7 mm in diameter, using the Energy button with Advanced Hemostasis. Do not attempt to seal vessels in excess of 7 mm in diameter.
- Blood and tissue buildup between the blade and shaft may result in abnormally high temperatures at the distal end of the shaft. To prevent burn injury, remove any visible tissue buildup at the distal end of the shaft.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eye wear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
- Do not attempt to bend, sharpen, or otherwise alter the shape of the blade. Doing so may cause blade failure and user or patient injury.
- To avoid user or patient injury in the event that accidental activation occurs, the instrument blade, clamp arm, and distal end of the shaft should not be in contact with the patient, drapes, or flammable materials while not in use.
- During and following activation in tissue, the instrument blade, clamp arm, and distal 7 cm of the shaft may be hot. Avoid unintended contact with tissue, drapes, surgical gowns, at all times.
- Do not introduce or withdraw the instrument with the jaws open through a trocar sleeve as this may damage the instrument.
- Care should be taken not to apply pressure between the instrument blade and tissue pad without having tissue between them. Clamping the tissue pad against the active blade without tissue on the full length of the blade will result in higher blade, clamp arm and distal shaft temperatures and can result in possible damage to the instrument. If this occurs, there may be an instrument failure, and the generator touchscreen displays a troubleshooting message.

- To avoid user or patient injury, do not activate an electrosurgical device in close proximity to the HARMONIC instruments. The aerosols created by the activation of the HARMONIC instruments in fatty tissue are potentially flammable.
- The entire exposed blade tip and any exposed blade shaft are active and will cut/coagulate tissue when the instrument blade is activated. Be careful to avoid inadvertent contact between all exposed blade surfaces and surrounding tissue when using the instrument.
- Use only the appropriate foot switch, instruments, and power cord to ensure that they are compatible with the generator.
- After removing the instrument, examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- Successful hemostasis may require adjunct measures when HARMONIC instruments are used to transect solid organs. Due to the difficulty of visualizing internal structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large vascular/biliary bundles when using the instrument under these conditions.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the HARMONIC system. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Incidental and prolonged activation against solid surfaces, such as bone, may result in blade heating and subsequent blade failure, and should be avoided.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization of single use devices may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.
- Reuse and improper reprocessing or resterilization of single use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

Disposal

Some internal components of the HARMONIC HD 1000i Shears instrument contain Lead (PZT Lead-Zirconate-Titanate). Disposal should be performed according to local requirements and regulations.

How Supplied

The HARMONIC HD 1000i Shears instrument is supplied sterile for single patient use. Discard the instrument after use.

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Ethicon Endo-Surgery, LLC ©Ethicon Endo-Surgery, Inc. 2016
Ethicon Endo-Surgery (Europe) GmbH
Ethicon Endo-Surgery, Inc.